

From: Toshihiko Kimura <t-kimura@asahi-spectra.co.jp>

Subject: Re: Thank you + a question

Date: 14 January 2025 at 4:45:54 GMT+2

To: Doron Chelouche <doron@sci.haifa.ac.il>

Dear Prof. Chelouche,

We made a simple video how to clean a surface of interference filter before.

<https://www.youtube.com/watch?v=9fVQJ-A08no>

We would like recommend that you prepare optical lens tissue and IPA (isopropyl alcohol).

And please blow any dusts on a surface completely.

Once again, I hope our narrow / medium band filters will work well at your telescope and be useful for your science!

Thank you,

Toshihiko